

Beakers, Berzelius, tall form



Code	Capacity (ml)	O.D. x H (mm)	Q'ty/Case
010040-5061A	50	44 x 75	6x10
010040-5071A	50	44 x 75	6
010040-10061A	100	50 x 90	6x10
010040-10071A	100	50 x 90	6
010040-20061A	200	60 x 111	6x6
010040-20071A	200	60 x 111	6
010040-30061A	300	65 x 135	6x4
010040-30071A	300	65 x 135	6
010040-50061A	500	77 x 150	6x4
010040-50071A	500	77 x 150	6
010040-100061A	1000	100 x 195	6x2
010040-100071A	1000	100 x 195	6
010040-200061A	2000	122 x 240	9
010040-300061A	3000	135 x 300	4

Material : High quality borosilicate glass (HARIO, Heat expansion coefficient $32 \times 10^{-7}/^{\circ}\text{C}$)