

Beakers, Griffin, low form with spout



Code	Capacity (ml)	O.D. x H (mm)	Q'ty/Case
010020-1071A	10	30 x 40	6
010020-2071A	20	34 x 45	6
010020-3071A	30	41 x 50	6
010020-5051A	50	46 x 61	6X12
010020-5071A	50	46 x 61	6
010020-10051A	100	55 x 73	6X6
010020-10071A	100	55 x 73	6
010020-20051A	200	67 x 90	6X6
010020-20071A	200	67 x 90	6
010020-30051A	300	78 x 105	6X4
010020-30071A	300	78 x 105	6
010020-50051A	500	92 x 121	6X2
010020-50071A	500	92 x 121	6
010020-100051A	1000	110 x 151	6X2
010020-100071A	1000	110 x 151	6
010020-200051A	2000	135 x 203	6
010020-300051A	3000	153 x 225	4
010020-500051A	5000	180 x 280	4
010020-1000051	10000	230 x 360	1

Material : High quality borosilicate glass (HARIO, Heat expansion coefficient $32 \times 10^{-7}/^{\circ}\text{C}$)