



Circulating Aspirators Model: WJ-15, WJ-20

The WJ-15 and WJ-20 are high-performance aspirators without connecting to a water faucet directly. They have much higher exhausting performance than traditional type. A unique construction not only prevents sucking the air but also suppresses noise with superior durability.

Features

- Stable exhaust
- High chemical resistance and durability
- Just turning on a switch to start operation without pump priming.



Specification

Code No.	044660-152	044660-202
Model	WJ-15	WJ-20
Performance	12 to 15L/min x 2pcs	
Motor	100/150W	
Aspirators	Metal water-jet x 2pcs	
Water bath	Polypropylene	
Suction port	2 places	
Water supply, Drain	Drain cock: nozzle ϕ 17mm Overflow: nozzle ϕ 17mm	Flexible hose: ϕ 12mm Overflow: nozzle ϕ 17mm
Accessories	G-14 hose nozzle 2pc, Junction tube(Y type) 1pc, Fuse(tube type) 1pc, Suction nozzle connector 2pcs	
Dimensions	255(W) x 360(D) x 415(H) mm	
Weight	Approx. 7.4kg	

* Possible to modify to a specified source voltage (modification or transforming).