



Dispenser, Manually type continuous "FINGER DISPENSER"



Code	Capacity(ml)	Range(ml)	Division(ml)	Accuracy(%)	Reproducibility(%)
025040-1	1	0.2 – 1	0.01	+/- 1	+/- 0.2
025040-2	2	0.4 – 2	0.01	+/- 1	+/- 0.2
025040-5	5	1 – 5	0.02	+/- 1	+/- 0.2
025040-10	10	2 – 10	0.1	+/- 1	+/- 0.2

※Tip is included 1 pc. Screw-cap bottle is not included.

Tips for FINGER DISPENSER

Code	Capacity(ml)	Quantity
025040-12A	1 – 2	5 pcs.
025040-510A	5 – 10	5 pcs.

Material: Polypropylene



SIBATA SCIENTIFIC TECHNOLOGY LTD.

1-1-62, Nakane Soka-City, Saitama, Japan

<http://www.sibata.co.jp/english/>

E-mail: overseas@sibata.co.jp

TEL:+81-48-933-1582 FAX:+81-48-933-1