

Economy Desiccators



EC-11



EUV-12



EUV-13



EC-21



EUV-22



EUV-23

Vertical Type

Code	047550-11	047550-12	047550-13
Model	EC-11	EUV-12	EUV-13
UV Protection	---	400nm	520nm
Shelf plate	3		
Inside dimensions	287(W) x 290(D) x 487(H)mm		
Outside dimensions	300(W) x 345(D) x 525(H)mm		

Horizontal Type

Code	047550-21	047550-22	047550-23
Model	EC-21	EUV-22	EUV-23
UV Protection	---	400nm	520nm
Shelf plate	2		
Inside dimensions	487(W) x 290(D) x 287(H)mm		
Outside dimensions	500(W) x 345(D) x 335(H)mm		