

Flasks, Erlenmeyer, narrow mouth



Code	Capacity (ml)	Neck I.D. mm	Max. O.D. mm	Hight mm	Q'ty/box
010530-1061A	10	13	30	50	10
010530-2061A	20	19	38	58	10
010530-3061A	30	19	46	65	10
010530-5051A	50	19	54	81	6 x12
010530-5071A	50	19	54	81	6
010530-10051A	100	20	66	102	6x8
010530-10071A	100	20	66	102	6
010530-20051A	200	27	80	126	6x4
010530-20071A	200	27	80	126	6
010530-30051A	300	31	90	148	6x4
010530-30071A	300	31	90	148	6
010530-50051A	500	37	108	171	6x2
010530-50071A	500	37	108	171	6
010530-100051A	1000	38	134	215	6
010530-200051A	2000	41	172	270	6
010530-300051A	3000	47	190	305	4
010530-500051A	5000	50	234	330	2

Material : High quality borosilicate glass (Heat expansion coefficient $32 \times 10^{-7}^{\circ}\text{C}$)