

Flasks, Round Bottom



Code	Capacity ml	Neck O.D. mm	Bulb O.D. mm	Overall H. mm	Q'ty/case
010500-5061A	50	20	53	105	36
010500-10041A	100	20	65	126	24
010500-20041A	200	25	80	155	24
010500-30041A	300	28	93	180	20
010500-50041A	500	28	108	210	20
010500-100011A	1000	34	135	260	6
010500-20001A	2000	38	170	325	6
010500-30001A	3000	40	200	370	4
010500-50001	5000	52	230	460	1
010500-100001	10000	70	290	580	1

Material : High quality borosilicate glass (HARIO, Heat expansion coefficient $32 \times 10^{-7}/^{\circ}\text{C}$)