

Flasks, filtering



Material : Borosilicate glass, JIS R3503 JR-2

Code	Capacity (ml)	Height (mm)	Lip dia. (mm)	Flask dia. (mm)	Inlet tube Dia.(mm)
017810-500	500	160	33	115	10
017810-1000	1000	195	36	142	10
017810-2000	2000	235	39	165	14
017810-3000	3000	255	39	192	14
017810-5000	5000	320	48	225	17

Material : High quality borosilicate glass, JIS R3503 JR-1

Code	Capacity (ml)	Height (mm)	Lip dia. (mm)	Flask dia. (mm)	Inlet tube Dia.(mm)
017800-500	500	160	34 ± 0.5	110	8
017800-1000	1000	204	35.5 ± 0.5	147	10