



Flasks, recovery , TS Joint



Code	Capacity ml	TS Joint	Bulb O.D. mm	Overall H. mm
005370-1510	10	15/25	32	65
005370-1525	25	15/25	42	80
005370-1550	50	15/25	49	99
005370-1950	50	19/38	50	101
005370-2450	50	24/40	50	103
005370-2950	50	29/42	50	105
005370-15100	100	15/25	64	124
005370-19100	100	19/38	63	126
005370-24100	100	24/40	63	128
005370-29100	100	29/42	63	130
005370-19200	200	19/38	76	156
005370-24200	200	24/40	76	158
005370-29200	200	29/42	76	151
005370-19300	300	19/38	85	176
005370-24300	300	24/40	85	178
005370-29300	300	29/42	85	180
005370-24500	500	24/40	100	203
005370-29500	500	29/42	100	205
005370-241	1000	24/40	136	232
005370-291	1000	29/42	136	228
005370-292	2000	29/42	167	280

Material : High quality borosilicate glass (HARIO, Heat expansion coefficient $32 \times 10^{-7}/^{\circ}\text{C}$)



SIBATA SCIENTIFIC TECHNOLOGY LTD.

1-1-62, Nakane Soka-City, Saitama, Japan

<http://www.sibata.co.jp/english/>

E-mail: overseas@sibata.co.jp

TEL:+81-48-933-1582 FAX:+81-48-933-1