



Graduated Cylinders, Custom A, with glass stopper (JIS R3505, Class A)



Specifications

Code No.	Capacity (ml)	Subdivision (ml)	Allowance (ml)	O.D. x H (mm)	Q'ty/case
026580-5A	5	0.1	±0.1	12 x 125	5
026580-10A	10	0.1	±0.2	13.1 x 190	5
026580-201A	20	0.2	±0.2	16 x 210	5
026580-251A	25	0.2	±0.25	17.4 x 220	5
026580-501A	50	0.5	±0.5	25 x 230	5
026580-1001A	100	1	±0.5	32 x 250	5
026580-2001A	200	2	±1.0	41 x 300	5
026580-2501A	250	2	±1.25	42 x 330	5
026580-5001A	500	5	±2.5	55 x 390	2
026580-10001A	1000	10	±5.0	70 x 450	2
026580-2000	2000	20	±10.0	90 x 530	1

※H(mm)=Without glass stopper

Material : HARIO, Heat expansion coefficient $32 \times 10^{-7}/^{\circ}\text{C}$



SIBATA SCIENTIFIC TECHNOLOGY LTD.

1-1-62, Nakane Soka-City, Saitama, Japan

<http://www.sibata.co.jp/english/>

E-mail: overseas@sibata.co.jp

TEL: +81-48-933-1582 FAX: +81-48-933-1