



## Magnetic Stirrer "THERMO MAGNESTIR"

### Model: MGH-110, MGH-320

The MGH-110 and MGH-320 can heat liquids up to 350°C with chemical-resistant ceramic plate.



Code No.	050620-110	050620-320
Model	MGH-110	MGH-320
Stirring capacity	Max. 1000ml in beaker	Max. 5000ml in beaker
Revolution speed	150 to 1500 rpm	
Speed control	Electronic control	
Motor	Shaded-pole motor, 1.8W	
Heater	0.4kW	1.1kW
Surface temp.	Max. 350°C	
Hot plate size	100(W) x 100(D)mm	180(W) x 180(D)mm
Magnet	Cobalt magnet	
Spin bar	PTFE coat φ8 x 35mm, 1pc	
Power supply	AC100V, 50/60Hz, 5A	AC100V, 50/60Hz, 12A
Dimensions	135(W) x 180(D) x 115(H) mm	220(W) x 245(D) x 115(H) mm
Weight	Approx. 1.8kg	Approx. 3.5kg

- It is possible to modify it to a specified source voltage (modification or transforming).