



Measuring Pipette, Custom A, Graduated to extreme tip (JIS R3505, Class A)



Specifications

Code No.	Capacity (ml)	Subdivision (ml)	Allowance (ml)	Length (mm)	Q'ty/case
021050-01A	0.1	0.001	±0.005	300	10
021050-02A	0.2	0.002	±0.005	300	10
021050-03A	0.3	0.002	±0.005	300	10
021050-05A	0.5	0.005	±0.005	300	10
021050-1A	1	0.01	±0.01	360	10
021050-2A	2	0.02	±0.015	360	10
021050-3A	3	0.02	±0.03	360	10
021050-5A	5	0.05	±0.03	366	10
021050-10A	10	0.1	±0.05	370	10
021050-20A	20	0.1	±0.1	435	10
021050-25A	25	0.1	±0.1	455	10

Material : High quality borosilicate glass (DURABAX, Heat expansion coefficient $49 \times 10^{-7}/^{\circ}\text{C}$)